


MANUAL OXYGEN RESPIRATOR

DIAGRAM OXYGEN-RESPIRATORYKIT



1. Oxygen cylinder 2 liter 200 bar.
2. Reducer 4,5 bar 50 liter/min.
Intergrated flow regulator
0-15 liter/min.
3. Resuscitator with inhalation
masks 3 & 5 (option)
4. Suction pump unit, manual
suction unit (option)
5. 2x disposable inhalation masks
with hose

CONTENTS:

- Suitcase
Dimensions: 50x40x20 cm
- Oxygen cylinder 2 liter 200 bar – Connection: G 5/8" Din
- Oxygen reducer 200 → 4,5 bar with manometer
- Guedel airway tube – size: 1
- Guedel airway tube – size: 2
- Guedel airway tube – size: 3
- Guedel airway tube – size: 4
- Disposable inhalation masks with hose and reservoir adult 2x
- Resuscitator with inhalation and oxygen bag (option)
- Connection hose for resuscitator
- Res-q-vac manual suction unit (option)
- Manual oxygen-respirator kit English

O2-Flow	Ventilation volume (ml) x frequency (cycly/min) oxygen flow (A)	
	Table 1	Table 2
A =	WITHOUT reservoir	WITH reservoir
Liter/min.	500 ml x 12	500 ml x 12
3	34%	68%
5	47%	85%
10	70%	100%

How to act in case victim is breathing:

Administer oxygen with inhalation mask(s)

1. Slowly open valve oxygen cylinder.
2. Set the oxygen flow regulator on the desired flow.
Adult: 6-8 liter/min. Child: 3-5 liter/min. Infant: 2-4 liter/min.
3. Remove viscous and solid food particles from mouth and troat by using finger, towel or suction unit. **Be aware of sweets chewing gum etc.**
4. Ease victim's head backwards to free air passage.
5. Place oxygen inhalation mask over the victims mouth and nose.
6. Gently press the mask to avoid oxygen leakage during inhalation.
7. After using equipment, close the valve of the oxygen cylinder and release the pressure, then close the flow regulator(s)
8. Replace used mask(s)
9. Refill cylinder.

How to act in case victim is NOT breathing:

Administer oxygen with resuscitator

1. Connect the oxygen tube between the flow regulator and resuscitator.
2. Connect eventual the oxygen reservoir to the resuscitator, see values in table 1 for use without reservoir and table 2 for use with reservoir.
3. Slowly open valve oxygen cylinder.
4. Set the oxygen flow regulator on the desired flow.
Adult: 6-8 liter/min.
Child: 3-5 liter/min.
Infant: 2-4 liter/min.
5. Remove viscous and solid food particles from mouth and troat by using finger, towel or suction unit. **Be aware of sweets chewing gum etc.**
6. Ease victim's head backwards to free air passage.
7. Place oxygen inhalation mask over the victims mouth and nose.
8. Gently press the mask to avoid oxygen leakage during inhalation.
9. Calmly squeeze into the balloon; Adults about: 12-15x per minute
Children/Infants about 20-25x per minute
10. After using equipment, close valve of the oxygen cylinder and release the pressure, then close the flow regulator(s)
11. Replace used mask(s) and balloon or clean balloon by sterilization.
12. Refill cylinder.

NOTE: flow settings are initial settings.

Administering oxygen should only be done authorized persons.